Mobile Docking System - M Series

User Guide

CE

SonoSite, Inc.

21919 30th Drive SE

Bothell, WA 98021

USA

T: 1-888-482-9449 or 1-425-951-1200

F: 1-425-951-1201

SonoSite Ltd

Alexander House

40A Wilbury Way

Hitchin

Herts SG4 0AP

UK

T: +44-1462-444800

F: +44-1462-444801

Caution: | United States federal law restricts this device to sale by or on the order of a physician.

SonoSite is a registered trademark of SonoSite, Inc.

Non-SonoSite product names may be registered trademarks or trademarks of their respective owners.

Patents pending.

P08326-01 11/2007

Copyright 2007 by SonoSite, Inc.

All rights reserved.

MDSm User Guide

Contents

Introduction	1
Safety	2
MDSm features	4
Connecting ultrasound system	6
Connecting printer and DVD recorder	6
Connecting printer only	9
Connecting external display	11
Connecting Triple Transducer Connect	13
Cables, connectors, and mini-dock symbols	14
Troubleshooting	19
Cleaning and disinfecting the MDSm	20
Cleaning the peripherals	20
Specifications	21

Introduction

The Mobile Docking System for M Series (MDSm) provides a platform for the M Series ultrasound system with connections for accessories and peripherals. This user guide provides information on the MDSm and on connecting accessories and peripherals.

For information on the peripherals, please see the manufacturer's operating instructions shipped with each product. For information on the ultrasound system, transducers, and additional safety information including warnings and cautions, see the ultrasound system user guide.

WARNING: To avoid the risk of patient injury, misdiagnosis, or user injury, read all the warnings in the ultrasound system user guide and supplemental guides.

Conventions, symbols, and terms

The user guide follows these conventions:

- A **WARNING** describes precautions necessary to prevent injury or loss of life.
- A **Caution** describes precautions necessary to protect the products.
- Numbered steps in procedures must be performed in a specific order.
- Bulleted lists present information in list format but do not imply a sequence.

Labeling symbols and terms used on the system are explained in the ultrasound system user guide.

Customer assistance and comments

Questions and comments are encouraged. SonoSite is interested in your feedback regarding the system and the user guide. Please call SonoSite at **1-888-482-9449**. For technical support, call **1-877-657-8118**, or if you are outside the USA, call the nearest SonoSite representative. You can also e-mail SonoSite at **comments@sonosite.com**.

Safety

For maximum safety, observe the following warnings and cautions before using the accessories and peripherals with the ultrasound system and transducers.

Electrical safety

WARNING:

To avoid the risk of electrical shock and fire hazard, a multiple portable socket outlet (MPSO) or extension cord should not be connected to the auxiliary mains socket on the MDSm.

To avoid the risk of electrical shock to the patient, do not touch the ungrounded connectors on the back of the mini-dock while simultaneously touching a patient.

To avoid the risk of electrical shock, use only properly grounded equipment. Shock hazards exist if the power supply is not properly grounded. Grounding reliability can only be achieved when equipment is connected to a receptacle marked "Hospital Only" or "Hospital Grade" or the equivalent. The grounding wire must not be removed or defeated.

To avoid the risk of electrical shock and fire hazard, inspect the power supply, plugs, cables, and power cords on a regular basis for damage.

To avoid the risk of electrical shock, use only accessories and peripherals recommended by SonoSite, including the power supply. Connection of accessories and peripherals not recommended by SonoSite could result in electrical shock. Contact SonoSite or your local representative for a list of accessories and peripherals available from or recommended by SonoSite.

To avoid the risk of electrical shock, use commercial grade peripherals recommended by SonoSite on battery power only. Do not connect these products to AC mains power when using the system to scan or diagnose a patient/subject. Contact SonoSite or your local representative for a list of the commercial grade peripherals available from or recommended by SonoSite.

WARNING: To avoid possible electrical shock or electromagnetic interference, verify proper operation and compliance with relevant safety standards for all equipment before clinical use. Connecting additional equipment to the ultrasound system constitutes configuring a medical system. SonoSite recommends verifying that the system, all combinations of equipment, and accessories connected to the ultrasound system comply with JACHO installation requirements and/or safety standards such as AAMI-ES1, NFPA 99 OR IEC Standard 60601-1-1 and electromagnetic compatibility standard IEC 60601-1-2 (Electromagnetic compatibility), and are certified according to IEC Standard 60950 (Information Technology Equipment (ITE)).

See the ultrasound system user guide for information and a declaration concerning Electromagnet Compatibility (EMC) for the ultrasound system.

Equipment safety

Caution: To avoid damaging the mini-dock, do not perform high current grounding impedance testing involving this part. Accessible metal of the mini-dock is not protectively earthed. For an explanation on the grounding of the mini-dock, see SonoSite Technical Bulletin 003/06 at www.sonosite.com. Select **Support and Service** and then select your system type.

Excessive bending or twisting of cables can cause a failure or intermittent operation.

Improper cleaning or disinfecting of any part of the accessories and peripherals can cause permanent damage. For cleaning and disinfecting instructions, see [“Cleaning and disinfecting the MDSm”](#) on page 20 and the manufacturer’s operating instructions.

Clinical safety

WARNING: To avoid misdiagnosis, do not use non-medical (commercial) grade peripheral monitors. These monitors have not been verified or validated by SonoSite as being suitable for diagnosis.

Electromagnetic Compatibility (EMC)

See the safety section of the ultrasound system user guide and supplemental guides for information on electromagnetic compatibility with IEC 60601-1-2, including the Manufacturer’s Declaration.

MDSm features

The MDSm provides a mobile work platform and provides a storage area for transducers and other supplies. The MDSm includes the mini-dock with connections to accessories and peripherals. Some cables are pre-routed on the MDSm. Other cables and cords must be installed manually. See ["Cables, connectors, and mini-dock symbols"](#) on page 14.

WARNING: To avoid the risk of electrical shock or injury, the equipment shall be connected to a center-tapped single phase supply circuit when users in the United States connect the equipment to a 240V supply system.

To avoid potential injury, do not touch the squared edges on the upper columns. The squared edges may be sharp.


Figure 1 MDSm: Front and Back Views

Table 1: MDSm System Features

Number	Description
1	Ultrasound system
2	Transducer, connector, and gel storage on both sides of MDSm
3	Side panels
4	Transducer storage
5	Height adjustment pedal
6	Locking wheels
7	Locking lever
8	Black-and-white printer bay
9	DVD recorder or color printer bay
10	Mini-dock connection panel (See “Cables, connectors, and mini-dock symbols” on page 14)
11	Cable cover
12	Upper columns
13	Storage
14	System power cord storage knobs
15	Back panel
<hr/>	
Adjust height	<ol style="list-style-type: none"> 1 Press down on the height adjustment pedal and grasp the side panels. 2 Pull up on the side panels to raise, and push down to lower.
Lock or unlock wheels	Press down on the locking levers on the front wheels to lock, and push up to unlock.

Insert system	1	Turn off the ultrasound system. The mini-dock may not function properly if you connect the system while the power is on.
	2	Raise the system locking lever.
	3	Slide the system into the MDSm.
	4	Lower the locking lever to secure the system.
<i>Note:</i> If at any time the ultrasound system should lose connectivity, turn off the system and reattach the mini-dock and system.		
Remove system	1	Turn off the ultrasound system.
	2	Close the lid on the ultrasound system.
	3	Lift up the locking lever and remove the ultrasound system.


Connecting ultrasound system

Connect AC power	1	Lock the wheels and remove the back panel.
	2	Ensure that the pre-routed DC power cord (A) is connected to the mini-dock and that the pre-installed AC power cord (I) is connected to the power strip.
	3	Connect the system AC power cord (K) to a hospital-grade electrical outlet.

Connecting printer and DVD recorder

WARNING:	<p>To avoid the risk of electrical shock or injury, an isolation transformer must be installed on the MDSm when the following conditions exist:</p> <ul style="list-style-type: none"> • More than two accessories or peripherals are used. • Both color and black-and-white printers are used. • Sony DVD recorder is used with any other accessory or peripheral. <p>Contact SonoSite or your SonoSite representative for specific information on the combinations of accessories and peripherals that require installation of an isolation transformer.</p> <p>To avoid potential injury, do not touch the squared edges on the upper columns. The squared edges may be sharp.</p>
Caution:	<p>To avoid damaging the system, use only peripherals recommended by SonoSite with the ultrasound system.</p>

The MDSm provides connections for a printer and a DVD recorder. See [Figure 2](#). For specific information on using the DVD recorders and printers, including warnings and cautions, see the manufacturer's operating instructions.

* Use adapter G2 and audio extension cables E1 with E.

** Use extension cable B1 and adapter B2 with B.

Figure 2 Connectivity of System, Printer, and a DVD Recorder

Connect cables	<ol style="list-style-type: none"> 1 Lock the wheels and remove the back panel. 2 Ensure that the pre-routed DC power cord (A) is connected to the mini-dock and that the pre-installed AC power cord (I) is connected to the power strip. 3 Connect the DVD recorder cables. See Figure 2 on page 7 and “Cables, connectors, and mini-dock symbols” on page 14. <ul style="list-style-type: none"> • AC power cord (I) • Pre-routed composite video cable (C) • Pre-routed S-video cable (B). Attach S-video adapter (B2) to extension cable (B1) and to S-video cable (B). • Pre-routed RS-232 cable (E). • Attach DVD adapter (G2) to the RS-232 connector, and attach audio extension cables (E1) to the red and white audio left/right input/output connectors. <p><i>Note:</i> RS-232 cable (E) is pre-routed through the MD5m column. You connect one end to the mini-dock. The other end has five connectors that you connect to the DVD recorder.</p> 4 Connect the printer cables. See Figure 2 on page 7 and “Cables, connectors, and mini-dock symbols” on page 14. <ul style="list-style-type: none"> • Pre-routed print control cable (H) • Printer AC power cord (I) • Composite video cable (D) 5 Connect system power cord (K) to a hospital-grade electrical outlet. 6 Install the back panel (vents at the top) by sliding the panel up and in.
Set up system	<ol style="list-style-type: none"> 1 Press the SETUP key on the ultrasound system. 2 Select Connectivity. 3 In the Printer list, select desired printer. 4 In the Video Mode list, select the appropriate video format: NTSC or PAL. 5 In the Serial Port list, select DVD. 6 Restart the system to activate DVD recorder connectivity.
Record images	<p><i>Note:</i> Do a test recording before a patient exam to ensure that all connectivity is correct.</p> <ol style="list-style-type: none"> 1 Turn on the DVD recorder. 2 Insert a DVD. 3 Do one of the following: <ul style="list-style-type: none"> • Use controls on the DVD recorder to start and stop recording. • If the A or B shortcut key is programmed to record, press it to start recording. Press it again to stop. <p>Each time you record, a new file saves to the DVD.</p> <p>By default the B shortcut key records. To program the A and B shortcut keys, see the ultrasound system user guide.</p>

- | | |
|-------------------------|---|
| Play back images | <ol style="list-style-type: none"> 1 Press the Function key (f) and then V on the system. 2 Press PLAY on the DVD recorder. Follow the manufacturer's instructions for playing back images. |
|-------------------------|---|

Connecting printer only

The MDSm provides connections for a black-and-white or color printer. For specific information on using the printer, including warnings and cautions, see the manufacturer's instructions.


Figure 3 Connectivity of System and a Printer on MDSm

- | | |
|-----------------------|---|
| Connect cables | <ol style="list-style-type: none"> 1 Lock the wheels and remove back panel.
If installing an older black-and-white or color printer, you may need to move padded screws to back holes. 2 Ensure that the pre-routed DC power cord (A) is connected to the mini-dock and that the pre-installed AC power cord (I) is connected to the power strip. 3 Connect the cables. See Figure 3 on page 9 and "Cables, connectors, and mini-dock symbols" on page 14. <ul style="list-style-type: none"> • Pre-routed print control cable (H) • Printer AC power cord (I) • Pre-routed composite video cable (C) 4 Connect system power cord (K) to a hospital-grade electrical outlet. 5 Install the back panel (vents at the top) by sliding panel up and in. |
| Set up system | <ol style="list-style-type: none"> 1 Press the SETUP key on the ultrasound system. 2 Select Connectivity. 3 In the Printer list, select the desired printer.
Only the printers listed are recommended for use with the system. |

Print screen

Note: Do a test printing before a patient exam to ensure that all connectivity is correct.

- 1 Ensure that the cables are correctly installed.
- 2 Turn on the printer.
- 3 (Color Printer Only) Change the input setting from S-video to video.
- 4 Do one of the following:
 - Use the controls on the printer for printing images. Follow the manufacturer's instructions.
 - If the A or B shortcut key is programmed to print, press it.
By default the A shortcut key prints. To program the A and B shortcut keys, see the ultrasound system user guide.

You can also save images and print them from the Patient List. For instructions, see the ultrasound system user guide.

Connecting external display

The MDSm provides connections for an external display. For specific information on using the external display, including warnings and cautions, see the manufacturer's operating instructions.


Figure 4 Connectivity of External Display on MDSm

- Connect cables**
- 1 Ensure that the pre-routed DC power supply cord (**A**) is connected to the mini-dock and that the pre-installed AC power cord (**I**) is connected to the power strip.
See [Figure 4](#) and [“Cables, connectors, and mini-dock symbols”](#) on page 14.
 - 2 Remove both protective covers from the DVI cable (**T**).
 - 3 Connect the following cables:
 - DVI cable (**T**)
 - AC power cord (**J**)
 - DC power cord (**S**)
 - 4 Feed the system AC power cord (**K**) through the power cord protector hole at back of the lower shelf, and then connect it to a hospital-grade electrical outlet.
A firm push may be required to fit the cord through the protector hole.

**Turn on
display and
adjust angle**

Adjust the display to the proper angle prior to use.

- 1 Turn on the display. (The switch is on the back.)
- 2 Turn on the ultrasound system.
- 3 Holding the display with both hands, rotate it forward or backward to the desired angle. Additional pressure may be required.


Figure 5 Adjust the Angle of the External Display

Caution:

A set screw controls the tightness of the display angle. To avoid damage to the display, do not adjust the set screw. Turning the set screw loosens the display, causing the display to drop.

Connecting Triple Transducer Connect

The Triple Transducer Connect (TTC) is an optional three-transducer module for the MDSm that allows connection of up to three transducers simultaneously while the ultrasound system is docked.

- | | |
|--|---|
| Connect transducer connector on TTC to system | <ol style="list-style-type: none"> 1 Raise locking lever. 2 Pull the transducer latch up and rotate 90°. 3 Align and insert the transducer connector to the connector on the bottom of the ultrasound system. 4 Turn the latch back 90° and press in, securing the transducer connector to the ultrasound system. 5 Place the system into the MDSm. <p>Connection is made when the locking lever is pushed down.</p> |
|--|---|


Figure 6 Connect Transducer Connector to System

- | | |
|----------------------------------|--|
| Connect transducer to TTC | <ol style="list-style-type: none"> 1 Pull the transducer latch up and rotate 90°. 2 Align and insert the transducer connector to the connector on the bottom of the TTC. 3 Turn the latch and press down, securing the transducer connector to the Triple Transducer Connect. |
|----------------------------------|--|

- | | |
|--------------------------|--|
| Select transducer | <p><i>Note:</i> When you use the TTC, there may be a minor reduction of transducer penetration.</p> <p>Push the transducer select button on the TTC to select the desired transducer. See Figure 6 on page 13.</p> <ul style="list-style-type: none"> • When selecting another transducer, the system restarts. The light on the button flashes momentarily and then turns solid green. • Only one transducer is accessible at a time. |
|--------------------------|--|

Cables, connectors, and mini-dock symbols

The cable kit includes cables for both the MDSm and the MDS Lite II. You need to use only the cables for your product.

Table 2: Connectivity Symbols on Mini-Dock

Symbol	Definition
	DC input
	Print control
	USB
	Ethernet
	RS-232 (DVD recorder or bar code scanner)
	S-video out
	S-video in
	DVI video out
	Composite video out
	Audio out
ECG	ECG/Footswitch

Table 3: Pre-routed Cables








Letter	Name	Illustration	Connector 1 [Illustration] [Connects to]	Connector 2 [Illustration] [Connects to]
A	DC power supply cord 6.8 ft./2 m		 Mini-dock	 Power supply
B	S-video cable (playback) 6 ft./1.8 m		 Mini-dock	 DVD recorder
C	Composite video cable 6 ft./1.8 m		 Mini-dock	 DVD recorder
E	RS-232 cable (control/audio) 6.5 ft./2 m		 Mini-dock	 DVD recorder
H	Print control cable 6 ft./1.8 m		 Mini-dock	 Printer

Table 4: Manually Installed Cables and Connectors

Letter	Name	Illustration	Connector 1 [Illustration] [Connects to]	Connector 2 [Illustration] [Connects to]	Product Used On
B	S-video cable (playback) 6 ft./1.8 m		 Mini-dock	 DVD recorder	MDS Lite II
B	S-video cable (record) 6 ft./1.8 m		 Mini-dock	 DVD recorder	MDS Lite II
B1	S-video extension cable 12 in./30.5 cm		 S-video adapter	 DVD recorder	MDSm
B2	S-video adapter		 S-video cable	 S-video extension cable	MDSm
C	Composite video cable 6 ft./1.8 m		 Mini-dock	 DVD recorder	MDS Lite II
D	Composite video cable 12 in./30.5 cm		 Printer	 DVD recorder	MDSm

Table 4: Manually Installed Cables and Connectors (Continued)






















Letter	Name	Illustration	Connector 1 [Illustration] [Connects to]	Connector 2 [Illustration] [Connects to]	Product Used On
E	RS-232 cable (control/audio) 6.5 ft./2 m		 Mini-dock	 DVD recorder	MDS Lite II
E1	Audio extension cables 12 in./30.5 cm		 RS-232 cable (audio output/input)	 DVD recorder	MDSm
E1	Audio extension cables 12 in./30.5 cm		 RS-232 cable (audio output/input)	 DVD recorder	MDSm
G2	DVD adapter 12 in./30.5 cm		 RS-232 cable	 DVD recorder	MDSm MDS Lite II
H	Print control cable 6 ft./1.8 m		 Mini-dock	 Printer	MDS Lite II
I	AC power cord 18 in./45.7 cm		 Power supply	 Power strip	MDSm MDS Lite II

Table 4: Manually Installed Cables and Connectors (Continued)

Letter	Name	Illustration	Connector 1 [Illustration] [Connects to]	Connector 2 [Illustration] [Connects to]	Product Used On
I	AC power cord 18 in./45.7 cm		 DVD recorder	 Power strip	MDSm MDS Lite II
I	AC power cord 18 in./45.7 cm		 Printer	 Power strip	MDSm MDS Lite II
J	AC power cord 39 in./1 m		 Printer	 Power strip	MDS Lite II
K	System AC power cord 10 ft./3 m		 Power strip	 Wall outlet	MDSm MDS Lite II (ships in separate box)
L	AC power cord 6 ft./1.8 m		 External display	 Power strip	MDSm
S	External display power supply		 Display	 Power supply	MDSm
T	External display DVI 3 ft./0.9 m		 Mini-dock	 Display	MDSm

Troubleshooting

Symptom	Solution
External display does not show an image.	<ul style="list-style-type: none"> Check the cable connections on display and stand. See Figure 4, "Connectivity of External Display on MDSm" on page 11. Ensure that the display is turned on and set up properly (i.e., ready to show correct video format: S-video, composite, or RGB). See the display manufacturer's instructions, if necessary.
External display does not project sound.	<ul style="list-style-type: none"> Check the cable connections on display and stand. See Figure 4, "Connectivity of External Display on MDSm" on page 11.
DVD recorder does not record.	<ul style="list-style-type: none"> Check the Serial Port setting in system setup. See "Set up system" on page 8. Check the cable connections on stand and DVD recorder. See Figure 2, "Connectivity of System, Printer, and a DVD Recorder" on page 7. Ensure that the DVD recorder is turned on and set up properly. See the manufacturer's instructions, if necessary.
DVD recorder does not play back.	<ul style="list-style-type: none"> Check the cable connections on display, stand, and DVD recorder. See Figure 2, "Connectivity of System, Printer, and a DVD Recorder" on page 7. Ensure that the display setting is set to S-video. Press the +/- key to toggle video format on the display (default setting is Composite).
Printer does not work.	<ul style="list-style-type: none"> Check the printer selection in system setup. See "Set up system" on page 9. Check the cable connections. See "Connecting printer and DVD recorder" on page 6 or "Connecting printer only" on page 9. Ensure that the printer is turned on and set up properly. See the printer manufacturer's instructions, if necessary.

Cleaning and disinfecting the MDSm

WARNING:	<p>To avoid electrical shock, remove the system and extra battery from the MDSm and unplug from power supply before cleaning.</p> <p>To avoid injury, always use protective eyewear and gloves when performing cleaning and disinfecting procedures.</p> <p>To avoid infection, if a pre-mixed disinfectant solution is used, observe the solution expiration date, and ensure that the date has not passed.</p> <p>To avoid infection, the level of disinfection required for a product is dictated by the type of tissue it contacts during use. Ensure that the solution strength and duration of contact are appropriate for the equipment. For information, see the disinfectant label instructions and the recommendations of the Association for Professionals in Infection Control and Epidemiology (APIC) and FDA.</p>
Caution:	<p>Use only recommended cleaners or disinfectants on surfaces. Immersion-type disinfectants are not approved for use on surfaces.</p>

Clean and disinfect	<p>The exterior surface of the MDSm can be cleaned and disinfected using a recommended cleaner or disinfectant.</p> <ol style="list-style-type: none">1 Remove the system, unplug the power supply, and detach any cables.2 Clean the surfaces using a soft cloth lightly dampened in a mild soap or detergent cleaning solution. Apply the solution to the cloth rather than the surface.3 Wipe the surfaces with the disinfecting solution. Theracide disinfectant is recommended.4 Air dry or towel dry with a clean cloth.
----------------------------	---

Cleaning the peripherals

For information about cleaning the peripherals, see the manufacturer's operating instructions.

Specifications

MDSm dimensions

Width: 22.6 in. (57.40 cm)
Length: 31.9 in. (81.03 cm)
Height: 34.2 - 46.2 in. (86.88–117.35 cm)
Height range: 12 in. (30.48 cm)
Weight: 78 lbs. (39.92 kg)

Temperature and humidity limits

Operating: 10–40°C; 15–95% R.H.
Shipping/Storage: -35–65°C; 15–95% R.H.

Peripherals

See the manufacturer's specifications for the peripherals.

Medical grade

Black-and-white printer
Recommended sources for printer paper: Contact Sony at 1-800-686-7669 or
www.sony.com/professional to order supplies or to obtain the name of the local distributor.
Color printer
DVD recorder

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